BANDE-CEINTURE ABDOMINALE

ABDOMINAL BAND BELT

6364

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EN | Instructions for Use

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SL | navodila za uporabo

SK | Návod na použitie

FR

BANDE-CEINTURE ABDOMINALE

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EN

ABDOMINAL BAND BELT

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

Without rigid back panel :

For abdominal conditions that may benefit from compression, such as:

- Abdominal strain or weakness of the abdominal area.

With rigid back panel :

For abdominal and low back conditions that may benefit from abdominal compression and gross motion restriction, such as:

- Abdominal strain or weakness of the abdominal area (post-partum, post-surgery) associated to : sciatica, lumbago / low back pain.

PERFORMANCE

The devices provide compression and / or support to the abdomen. The devices reduce low back pain.

The devices provide lumbar spine gross motion restriction, when the removable back panel is applied.

CONTRAINDICATION

Must not be used by individuals for whom compression is contraindicated.

MEASUREMENT AND SIZE CHART

The size is to choose in accordance with the waist circumference (in cm), measured 2 cm above the iliac crests.

	1	2	3	4	5
	78-88	88-98	98-108	108-118	118-128

WARNINGS • PRECAUTIONS OF USE

Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.

This product contains Naturel Rubber Latex which may cause allergic reactions.

Do not wear for extended periods without medical follow-up. The orthosis must not be used in direct contact with a damaged skin. The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional. The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Gradual withdrawal possibly combined with functional physiotherapy.


FEATURES OF THE DEVICE

This device is a band-belt stretchable in the lengthwise direction. The band-belt has a standard height of 25 cm and covers the abdomen. The belt is adjustable by adjusting the length of the belt with a hook-and-loop fastening system. An adjustable hypogastric strap is placed on the front of the band-belt. The belt has also a removable back panel composed of four stays to support the lumbar region.

EN

ABDOMINAL BAND BELT

EN

SET UP AND OPERATION

The first application should be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

If the back panel is not needed:

1. Remove the back panel.
2. Unfasten the two tabs **1D** & **1G** of hypogastric strap from their velcro fastenings.
3. Position the seam **2** of the tabs over the belly button.
4. Adjust and attach the velcro on the side of the belt to your waist **3**, making sure to keep the seam **2** of the two tabs at the belly button.
5. Then attach the 2 tabs **1D** & **1G** of hypogastric strap to the desired tightness, beginning with the left tab.

If the back panel is required:
6. Follow steps 2 and 3.
7. Make sure that the back panel is centered on the spine, so that it is in line with the seam **2**.
8. Follow steps 4 and 5.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Hand wash in warm soapy water (maximum 40 ° C), Rinse thoroughly. Wring out gently; Do not bleach; Dry flat, do not tumble dry; Do not iron; Do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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