GIBAUD

COUDIÈRE

ELBOW SLEEVE

6728

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EN | Instructions for Use

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

For elbow conditions that may benefit from compression, such as:

- Medial and lateral epicondylitis.
- Elbow pain and swelling.
- Elbow instabilities.

PERFORMANCE

The devices provide compression to the elbow.

CONTRAINDICATION

The brace must not be used by individuals for whom compression is contraindicated.

WARNINGS • PRECAUTIONS OF USE

This product contains Naturel Rubber Latex which may cause allergic reactions.

The device should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

To fit, pass the hand through the elbow sleeve.

Once in position, the logo should be at the end nearest to the hand on the forearm (not at the end nearest to the shoulder).

The openwork knitted section must be in the fold of the elbow.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Store the sleeve away from sources of extreme heat (including radiators and tumble dryers).

Hand wash in warm soapy water (maximum 40°C), rinse thoroughly; wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

The size is to choose in accordance with the circumference of the forearm below the elbow (in cm)

| | | | | | |
|--|-------|-------|-------|-------|-------|
|  | 1 | 2 | 3 | 4 | 5 |
| | 22-24 | 24-26 | 26-28 | 28-30 | 30-32 |

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GIBAUD S.A.S.U