GIBAUD

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
If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopaedic consultant for advice.

INDICATIONS

For knee conditions that may benefit from compression, such as :

- Patella tendinopathy.
- Patellofemoral Osteoarthritis (PFOA).
- Chondromalacia.
- Mild sprains and strains.
- Knee pain.
- Feeling of instability.

PERFORMANCE

The device provides realignment of the patella.

The device provides class 3 compression.

The device provides proprioception.

CONTRAINDICATION

Must not be used by individuals from whom compression is contraindicated.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

Users with medical conditions that decrease blood circulation in limbs (e.g. Diabetes or peripheral vascular diseases) should consult their doctor before using the brace.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Elastomeric patellar ring with groove for the patellar tendon.

Steel spiral stays.

SET UP AND OPERATION:

The first application should be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

When fitting, the kneecap should be positioned in the centre of the kneecap ring, and care should be taken to position the stabilization groove over the tibial tuberosity.

CARE INSTRUCTION AND MAINTENANCE

Hand wash in warm and soapy water (maximum of 40°C); rinse thoroughly; wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the circumference of your leg measured (in cm) around the kneecap and 15 cm above especially for sizes 4, 4+, 5 et 5+, with the knee slightly flexed.

1	2	3	4	4+	5	5+
36-41	41-46	46-51	51-56	56-61	56-61	61-66
28-33	33-38	38-43	43-48	43-48	48-53	48-53

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GIBAUD S.A.S.U