MALLÉOGIB® 3D

6679

FR

MALLÉOGIB® 3D

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EN

MALLÉOGIB® 3D

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopaedic consultant for advice.

INDICATIONS

For Ankle conditions that may benefit from compression, such as;

- Ankle injury rehabilitation
- Mild sprains and strains
- Mild instability

PERFORMANCE

- The device provides ankle support
- The device provides proprioception
- The device provides class 3 compression

CONTRAINDICATION

Must not be used by individuals from whom compression is contraindicated

WARNING • PRECAUTION OF USE

- The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.
- Do not wear for extended periods without medical follow-up.
- The orthosis must not be used in direct contact with damaged skin.
- If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.
- The orthosis is a single patient use device. Do not use for multiple patients
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Competent Authority of the Member State is identified on this web site: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

FEATURES AND BENEFITS

- Latex free

SET UP AND OPERATION

- In addition to taking measurements, testing is required for a proper fitting, and for the patient to learn the proper placement.
- The ankle sleeve can be used on the left or right ankle.
- The ankle sleeve is put on by sliding it over the ankle like a sock. The Gibaud logo should be on the top of the ankle. -Place the lower edge (without logo) under the arch of your foot.

CARE INSTRUCTION AND MAINTENANCE

Hand wash in warm and soapy water (maximum of 40°C); rinse thoroughly; wring out gently; Do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean. The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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
EN | Instructions for Use

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MEASUREMENT AND SIZE CHART

The circumference of your ankle should be measured above the malleoli.

1	2	3	4	5
18-20	20,5-23	23,5-26	26,5-29	29,5-32

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