GIBAU

MANUGIB® TENDINITE

MANUGIB® TENDONITIS

6686

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EN | Instructions for Use

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
MANUGIB® TENDONITIS

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

For conditions that may benefit from restricting motion of the wrist, such as:

- Wrist tendinopathy.
- Mild to moderate sprains and strains.
- Post-cast healing.
- Bursitis.
- Arthritis.
- Other selected soft tissue injuries.

PERFORMANCE

The brace shall provide motion restriction of the wrist.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

When adjusting the strap, tighten it so that it is firm but comfortable.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The orthosis is a single patient use device. Do not use for multiple patients.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Single adaptation tab.

SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

During fitting, the applicator must shape, if necessary, the splints of the orthosis to adapt them to your morphology.

1. Open the main flap.
2. Put the brace on like a glove;
3. Center the palm splint in the hand;
4. Tighten with the flap to obtain a firm but comfortable hold.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the device away from sources of excessive heat.

Hand wash in warm soapy water (maximum 40°C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the circumference of the wrist, below the styloid (in cm).

	1G - Left / 1D - Right	2G - Left / 2D - Right	3G - Left / 3D - Right
	13,5-15,5	15,5-18	18-21

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GIBAUD S.A.S.U