GIBAU

MANUGIB® TRAUMA POINET-POUCE

MANUGIB® TRAUMA WRIST AND THUMB

6650

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EN | Instructions for Use

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
MANUGIB® TRAUMA WRIST AND THUMB

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

For conditions that may benefit from restricting motion of the wrist and thumb, such as:

- Wrist/thumb trauma.
- Post-cast healing.
- Post-operative use.
- De Quervain's tenosynovitis.
- Gamekeeper's thumb.
- Wrist/thumb overuse injury.
- Other selected soft tissue injuries.

PERFORMANCE

The brace shall provide motion restriction of the wrist and thumb.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

When adjusting the strap, tighten it so that it is firm but comfortable.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Padd of comfort fills the palm arch.

4 tightening straps.

Large palm splint and thumb splint.

Removable conformable splints.

SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

During fitting, the applicator must shape, if necessary, the splints of the orthosis to adapt them to your morphology.

- 1) Loosen the straps.
- 2) Slip the brace on like a glove, positioning the stay on the palm side of the hand and putting the thumb into the space provided.
- 3) Tighten the wrist straps one after the other until the fit is snug but comfortable; the thumb splint must be positioned without any painful pressure points.
- 4) Then tighten the thumb strap to obtain a firm but comfortable hold.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the orthosis away from sources of excessive heat.

Hand wash in warm soapy water (maximum 40°C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the circumference of the wrist, below the styloid (in cm).

	0G - Left / 0D - Right	1G - Left / 1D - Right	2G - Left / 2D - Right
	10-14	18-22	18-22

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GIBAUD S.A.S.U