GIBAUD

ATTELLE DE DOIGT POIGNET FINGER AND WRIST SPLINT

6610

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EN | Instructions for Use

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

Trauma to the long fingers that may benefit from immobilization, such as:

- Sprains.
- Strains.
- Dislocation.
- Stable fractures.

PERFORMANCE

The brace shall provide immobilization of the long fingers.

The brace shall provide effective treatment for long finger injuries.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires absolute immobilization or surgical procedure. The brace should be fitted by a healthcare professional.

When adjusting the strap, tighten it so that it is firm but comfortable.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional. The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Adaptable splint according to medical recommendation.

Possible positioning in dorsal or palmar depending on the location of the trauma.

Immobilization in syndactyly: 2nd with 3rd finger and 4th with 3rd finger.

SET UP AND OPERATION

The first application should be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

Depending on the injury and medical advice, the splint may be fitted in the palm position (1) or back-of-the-hand position (2-3).

Immobilization position:

In the absence of medical advice to the contrary, immobilization should be in the protective position (immobilization in the intrinsic plus position) to prevent stiffening (4):

- Metacarpophalangeal joint (MCP): 60° - 70°.
- Proximal interphalangeal joint (PIP): 10° - 20°.
- Distal interphalangeal joint (DIP): extension position.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the orthosis away from sources of excessive heat.

Hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the length of the hand, measured between the wrist and the end of the injured finger (in cm).

1	2	3	4
< 13	13 - 17	17 - 20	20 - 23

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GIBAUD S.A.S.U