## ACTION V

FR

## ACTION V

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EN

## ACTION V

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

## INDICATIONS

Lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain.
- Lumbosacral strains.
- Chronic low back pain.
- Sciatica.

## PERFORMANCE

The device provides lumbar spine gross motion restriction. The device reduces stress on the low back in flexion and extension

## CONTRAINDICATION

Must not be used by individuals for whom compression is contraindicated.

**WARNING • PRECAUTION OF USE**  
**Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.**

The orthosis must not be used in direct contact with a damaged skin. This product contains Naturel Rubber Latex which may cause allergic reactions.


Do not wear for extended periods without medical follow-up. The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional. The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Gradual withdrawal possibly combined with functional physiotherapy.

## MEASUREMENT AND SIZE CHART

The device is available in two heights: 21 cm (6673) & 26 cm (6672). The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm)

## FEATURES AND BENEFITS

The device is a lumbar stretchable belt featuring four V-shaped posterior stays.  
Ergonomic finger loops facilitate installation.

## SET UP AND OPERATION

The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user, a family member or a friend. During testing, the fitter must shape the V-shaped ribs of the brace to adjust it to your form.  
- Wear over an undergarment or light clothing.

1) Position the bottom part of the belt at the top of the gluteal cleft  
2) Center the dorsal stays on either side of the spine. The Gibaud logo must be on the right side; the stays form a V on either side of the lumbar column.  
3) The belt must be closed by stretching the two sides symmetrically using the finger holes.

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