**GIBAUD**

# **GENUGIB® POST-OP**

**6498**

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**EN** | Instructions for Use

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**If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.**

## **INDICATIONS**

Post-op/post-injury knee immobilization and / or Range Of Motion (ROM) control.

## **PERFORMANCE**

The device provides motion restriction of the knee joint.

The device provides variable settings of ROM in the sagittal plane.

## **CONTRAINDICATIONS**

None known.

## **WARNINGS • PRECAUTIONS OF USE**

The brace must not be used in cases where the injury requires absolute immobilization or surgical procedure. When adjusting the strap, tighten it so that it is firm but comfortable.

Range of motion adjustment, as prescribed, should be done by a qualified person before application of the brace.

If any problems occur with the application of the orthosis, contact your healthcare professional.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and

contact your healthcare professional. The orthosis is a single patient use device. Do not use for multiple patients.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **SET UP AND OPERATION**

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

During fitting, the applicator must, if necessary, adapt the hinges of the orthosis to your morphology.

Adjustment of the ROM is absolutely essential before placement. The brace is delivered blocked in 15° flexion.

### **1) For adjusting the ROM:**

To change the setting, unscrew the transparent cover with a small coin, and place the 4 pins in the orifices facing the prescribed angles.

## **FEATURES OF THE DEVICE**

Polycentric hinges.

Anterior Tibial Tuberosity (ATT) silicone insert.

Amplitude adjustment system.

FLEXION	0°	15°	30°	50°	70°	90°	110°
EXTENSION	0°	15°	30°	50°	70°	90°	110°

The orthosis can be worn over trousers.

**2) For the setting of the orthosis:**

1. Open the upper and lower parts of the brace.
2. Pull the brace on, with the label positioned above, so the hinge is at knee level and the silicone pad is below the patella on the anterior tibial tubercle when the lower straps are closed. The two side cushions are detachable and can be removed in the case of discomfort.
3. Finish adjustment by positioning the sidebars to obtain firm but comfortable support. The hinges must be parallel and centered on the middle of the knee.

### CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Store the brace away from sources of excessive heat.

Remove the lateral hinges then hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from sources of excessive heat (ex: radiator), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### MEASUREMENT AND SIZE CHART

The brace is in a single size. The product is adjustable in width.

Height (cm)	50
Size	1

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GIBAUD S.A.S.U