GIBAUD

GENUGIB[®] ROM

6621

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EN | Instructions for Use

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

For knee conditions that may benefit from increased stability around the knee, such as:

- Mild to moderate ACL or PCL, MCL and LCL instabilities, sprains and strains.

PERFORMANCE

The device shall provide increased stability around the knee.

The device provides limitation in flexion and extension of the knee joint.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization by cast or resin or surgical procedure.

Do not wear for extended periods without medical follow-up.

When adjusting the straps, tighten them so that it is firm but comfortable.

The orthosis must not be used in direct contact with a damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional. The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Polycentric hinges.

Hinge locking and adjustment according to medical recommendations.

Anatomical shape achieved by tightening six opposing straps above and below knee.


Opens completely from the front.

Patellar opening with semi-circular rings and foam pad to hold the anterior tibial tubercle in place.

Removable condyle pads that can be adjusted to fit the morphology.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the circumference of the knee (in slightly bent position), measured at the center of the patella (in cm).

1	2	3	4	5	6
31-34	34-37	37-40	40-43	43-46	46-50

SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

1) For adjusting the range of motion: It must imperatively be carried out by the healthcare professional before fitting the orthosis.

To adjust, slide the hinge cap downwards and remove it.

Position the red angle stops (extension and flexion) depending on the range of motion or blocking recommended by the prescriber. **Example of adjustment using the red stops.**

Close the cap again by sliding upwards.

2) For the setting of the orthosis: Open the brace up completely at the front.

Place the two grey patellar half-rings on either side of the patella.

Fasten the bottom then the top using the self-fastening tabs.

Position the silicone pad below the kneecap on the tibial crest and fasten the straps.

Tighten the straps until the fit is firm but comfortable.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep brace away from direct heat sources.

Hand wash after ribs have been removed in warm soapy water (maximum 40° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from direct heat sources (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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GIBAUD S.A.S.U