GIBAUD

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

For conditions that may benefit from added support to the ankle, such as:

- Ankle injury rehabilitation/ immobilization (Grade I/II sprains).
- Support chronic instability.
- Prophylactic use.

PERFORMANCE

The device provides functional support.

The device increases proprioception. The device prevents (further) injuries to the ankle.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

This product is designed to be worn over a sock.

Inspect your foot, ankle and this product before and after each use.

Do not over-tighten the brace, risking loss of blood flow. Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease need to be extra careful and are advised to consult with their healthcare professional before using the brace.

The brace must not be used in direct contact with a damaged skin.

Never wear the brace directly over an open wound.

This brace has been designed and tested based on single patient usage and is not recommended for multiple patient use.

If any problems occur with the use of the brace, e.g. pain or appearance of local signs immediately contact your healthcare professional.

If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white or cold) while using this brace, discontinue the use of it immediately and consult with your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Manual tightening of lace having it firm but comfortable.

Rigid morpho-adaptable parts.

SET UP AND OPERATION

The first application should be carried out by a healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend. For optimal comfort, the brace should be worn over a sock.

1) Insert the foot in the brace, with the heel in the heel opening.

2) Tighten the laces of the first four holes so that they are firm but comfortable.

3) Complete the fitting with the quick-tightening system, as you would when lacing up a shoe.

4) Adjust the tongue if necessary.

5) Slide the fitted foot into a spacious, stable shoe: the tightening of the shoe should be adjustable.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the ankle circumference (in cm).

1	2	3
19-22	23-26	27-29

6) Adjust the tightness of the shoe. Loosening may occur with use, re-tighten if necessary. If any localised discomfort or pinching occurs, particularly when used for sporting activities, loosen the laces and reposition the brace.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

After having removed the laces, hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a heat source (ex: radiator), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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