SOUTIEN ROTULIEN PATELLAR SUPPORT

5004

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

- Patellar tendinitis / jumper's knee.
- Instability of the patella.

PERFORMANCE

The device provides proprioception. The device provides realignment of the patella. The device provides stabilization of the patella. The device provides pain reduction.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

When adjusting the strap, tighten it so that it is firm but comfortable.

If any problems occur with the application of the orthosis, contact your healthcare professional.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact healthcare professionals.

The orthosis is a single patient use device. Do not use for multiple patients.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

A double central pad. A closing strap.

SET UP AND OPERATION


The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

Place the double central pad at the level of the patellar tendon and adjust the strap so as to obtain a firm but comfortable tightening.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Hand wash in warm soapy water (maximum 40 °C), rinse thoroughly. Wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the circumference of the knee under the patella (in cm).

Size	0	1
Knee circumference under the kneecap (en cm)	< 31	> 31

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