**GIBAUD**

# EQUALIZER® AIR

6715

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**EN** | Instructions for Use

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

## **INDICATIONS**

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains.
- Stable fractures.
- Soft tissue injuries.
- Post-operative, trauma or rehab use.

## **PERFORMANCE**

The device provides immobilization/motion restriction.

The device has rocker bottom that provides natural and stable gait.

The device has shock absorbing sole.

## **CONTRAINDICATION**

No absolute contraindications.

## **WARNINGS • PRECAUTIONS OF USE**

Do not remove the device unless instructed by physician.

This product has been designed and tested based on single patient usage and is not recommended for multiple patient usage.

The orthosis must not be used in direct contact with a damaged skin.

If skin irritations of any kind develop, consult healthcare professional immediately.

Be sure the patient performs circulation checks. If a loss of circulation is felt, or if the device feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.


When prescribed to patients with diabetes, vascular deficiencies and neuropathies extra caution and attention to patient is needed.

Do not remove the inner foam pads which are adhered to the plastic walker boot.

This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your healthcare professional.

## **FEATURES AND BENEFITS**

- A. Pull tabs
- B. Adjustable air pump
- C. Walker struts
- D. Adjustable pneumatic softgood
- E. Wide rocker bottom base
- F. Adjustable contact closure straps
- G. Foam pads

## **MEASUREMENT AND SIZE CHART**

The shoe size is to choose in accordance with the length of your foot.

Size	0	1	2	3	4
Shoe sizes (UK)	2-3.5	3.5-6	6-9.5	9.5-11	11-13

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

1. Remove softgood by opening contact closure straps and lifting out of brace. Do not remove pull tabs at this time.
2. Place lower leg and foot into softgood. Close softgood snugly over leg.
3. Insert foot into brace by gently spreading uprights. Do not remove pull tabs at this time.
4. Ensure struts are centered over the ankle bone (malleoli). The lower leg should be kept at a 90° angle with the foot.
5. Remove pull tabs located in between struts and softgood. This allows the softgood to be secured into position.
6. To accommodate for rubbing or spacing, additional foam pads may be placed inside the softgood, in between the softgood and plastic, or under straps for additional comfort.
7. Secure contact closure straps, starting at the toes and working upward. For best fit, secure straps in the correct order. Strap adjustments may be necessary as foot, ankle and leg swelling patterns change.
8. Inflate air chamber by turning knob clockwise until secure and pressing and releasing the air bulb repeatedly. Continue until softgood feels stable and is cushioning the lower leg. Do not over inflate.
9. Air can be released from air chamber by turning the inflation/deflation knob counterclockwise. To re-inflate, repeat step 8.

10. When prescribed by physician, remove leg from product by releasing contact closure straps and opening softgood.

Pull foot and leg forward and out, leaving softgood inside the brace. If softgood is removed, follow Steps 1 through 7 to reapply.

## CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Remove softgood from walker plastic before washing.

Hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a heat source (ex: radiator), do not tumble dry; do not iron; do not dry clean.

Walker plastic can be washed with a damp cloth and mild soap.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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