G I B A U D

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**If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.**

### INDICATION:

For low back conditions that may benefit from gross motion restriction, such as;

- Sciatica, Cruralgia in acute phase
- Lumbar traumatism (mild settlement fracture)
- Debilitating chronic low back pain
- Post-surgical stabilization

### PERFORMANCE

The devices provide lumbar spine gross motion restriction

The devices reduce low back pain

### CONTRAINDICATION

Must not be used by individuals for whom compression is contraindicated.

### WARNING • PRECAUTION OF USE

Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.

The orthosis must not be used in direct contact with a damaged skin

This product contains Naturel Rubber Latex which may cause allergic reactions.

Do not wear for extended periods without medical follow-up

The orthosis is a single patient use device. Do not use for multiple patients.

This product is sensitive to heat: never expose it to sunlight, place it on a radiator or attempt to make alterations without the help of the applicator

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional  
The user and/or patient must report any serious incident that has

occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Gradual withdrawal possibly combined with functional physiotherapy.

### FEATURES AND BENEFITS

This device is made of a 3.2 mm-thickness thermoformable perforated plastic back panel (posterior rigid plastic shell) up to D6.

The anterior plastron (front pad) has adjustable tightening straps and there are two shoulder straps exerting extension bracing. Under-thigh straps ensure apparatus stability.

### SET UP AND OPERATION

A fitting is compulsory for a correct adaptation and a good understanding of the installation.

The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

### Adaptation:


The adaptation must be carried out by a qualified applicator.

Possible retouching needs can then be highlighted: support at the floating ribs, hyper support or detachment at the sacrum ...

It is possible to adapt the plastic by modifying the curvature after heating with a hair dryer or a heat gun (sweeping between 1 to 3 minutes depending on the power of the device. It is not recommended to modify the central part of the corset (stiffening curve).

### MEASUREMENT AND SIZE CHART

The size is to choose in accordance with the waist circumference ), measured 2cm above the iliac crests (in cm)

1	2	3
70-85 cm	85-100 cm	100-115 cm

**Set up:**

This device is intended to be worn between clothing and underwear.

The plastic frame must be worn over a cotton undergarment.

1) Unhook the under-thighs from the plastron

- open the front plate by the central closure

- put on the shoulder straps

2) Close the central part of the faceplate by distributing moderate pressure

3-4) Adjust and distribute the compression using the 3 tightening straps (2 on the left and 1 on the right)


5) Adjust the recall of the suspenders thanks to their hooked ends on the faceplate

6) Position the under-thighs by fixing the hook fabric part on the velvet fabric of the plastron

1

2

**CARE INSTRUCTION, MAINTENANCE AND DISPOSAL**

Store the brace away from sources of extreme heat (including radiators and tumble dryers). Hand wash in warm soapy water (maximum 40 ° C); rinse thoroughly. Wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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GIBAUD S.A.S.U