G I B A U D

# GENOUILLÈRE

# KNEE SLEEVE

6726, 6733

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**EN** | Instructions for Use

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**If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopaedic consultant for advice.**

## INDICATIONS

For knee conditions that may benefit from compression, such as;

- Feeling of instability.
- Mild sprains and strains.
- Knee pain.

## PERFORMANCE

The device provides class 3 compression. The device provides proprioception.

## CONTRAINDICATION

Must not be used by individuals from whom compression is contraindicated.

## WARNINGS • PRECAUTIONS OF USE

This product contains Naturel Rubber Latex which may cause allergic reactions. The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

Users with medical conditions that decrease blood circulation in limbs (e.g. Diabetes or peripheral vascular diseases) should consult their doctor before using the brace.

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

When fitting, place the patella positioning oval over the kneecap.

## CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Hand wash in warm and soapy water (maximum of 40°C); rinse thoroughly; wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## MEASUREMENT AND SIZE CHART

The size of the orthosis is to be chosen according to the circumference of the knee (in slightly bent position), measured at the center of the patella (in cm).

1	2	3	4	5
28-33	33-38	38-43	43-48	48-53

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GIBAUD S.A.S.U