GIBAUD

# LOMBOGIB® UNDERWEAR

6532/6655/6664/6663

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**EN** | Instructions for Use

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**SV** | Bruksanvisning

**NO** | Bruksanvisning

**DA** | Brugsanvisning

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

## INDICATIONS

Lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain.
- Sciatica.
- Lumbosacral strains.
- Chronic low back pain.

## PERFORMANCE

The device provides lumbar spine gross motion restriction.

The device reduces stress on the low back in flexion and extension

## CONTRAINDICATION

Must not be used by individuals for whom compression is contraindicated.

## WARNING • PRECAUTION OF USE

Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.

The orthosis must not be used in direct contact with a damaged skin.

Do not wear for extended periods without medical follow-up.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Gradual withdrawal possibly combined with functional physiotherapy.

## FEATURES AND BENEFITS

The device is a lumbar stretchable belt featuring four V-shaped posterior stays.

Ergonomic finger loops facilitate installation.

## SET UP AND OPERATION

The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user, a family member or a friend.

- The belt is worn under the clothes or over the undergarments.

1) Position the bottom of the belt on the top of gluteal cleft.

2) Centre the dorsal rib on either side of the spinal column. The Gibaud logo must be on the right.

3) Tighten the belt based on the intensity of the pain and type of activity, to obtain a tighten firm but comfortable.

## CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

- Keep the belt away from direct heat sources (radiators, dryers...).

30°C mild fine wash (Hand wash or machine wash, mild process); rinse thoroughly. Wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations

## MEASUREMENT AND SIZE CHART

The device is available in two heights: 21 cm (6662/6663) & 26 cm (6532/6655)  
The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm)

	0	1	2	3	4	5
	70-80 cm	80-90 cm	90-100 cm	100-110 cm	110-125 cm	125-140 cm

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GIBAUD S.A.S.U